

## Item data sheet



### Butter-Croissant with apricot filling

<b>Item no.</b>	<b>922</b>
Designation of the food item	Croissant with apricot filling
EAN	4009837009220
Convenience level	preproved
This item is vegetarian.	

### Units

Unit	PIECE *	BAG	CARTON	LAYER	PALLET
Quantity per base unit	1	30	60	480	3.840
Gross weight <sup>1</sup>	0,09813	2,9439	5,8878	47,1024	376,8192
Net weight	0,09	2,70	5,40	43,20	345,60
L x W x H (mm)	150 x 65 x 40	520 x 350 x 0	396 x 296 x 242	-	1.200 x 800 x 2.086
Diameter (mm)	-	-	-	-	-
Cartons/pallets	-	-	-	-	64

\* Base unit of the item

<sup>1</sup> Weight in frozen state incl. Packaging

### Measurements

Measurements: l 15 x w 6.5 x h 4 cm

### Nutritional values

### Ingredients

### Allergens

#### Butter-Croissant with apricot filling

energy	1330 kJ / 317 kcal	WHEAT FLOUR, water, BUTTER (15%), sugar, glucose syrup, apricot pieces (4%), yeast, WHEAT GLUTEN, apples, salt, shea nut fat, SKIMMED MILK, rapeseed oil, dextrose, WHEAT STARCH, vegetable protein, MILK PROTEIN, colouring agent lutein, preservatives potassium sorbate, antioxidizing agent ascorbid acid, acidulant regulator(citric acid, sodium citrate), thickening agent pectin, emulsifier	Contains: GLUTEN AND - PRODUCTS, MILK AND - PRODUCTS
fat	13,5		
thereof saturated fatty acids	8,4 g		
carbohydrates	42,3 g		
thereof sugar	16,0 g		
protein	6,1 g		

salt	0,87 g	ditartaric acid ester from mono and di-glyceride of fatty acids, aroma, enzyme May contain: EGG AND - PRODUCTS, SOYA, NUT PRODUCTS	
------	--------	---	--

All data refer to 100 g

### Baking instructions

30 minutes defrost and bake for 20-24 minutes at 180-220°C with steam.

### Best before date (from date of production)

12 months

### Transport and storage conditions

This product must be transported and stored at a temperature of at least -18°C.

### Gene Ordinance

This product does not contain any genetically modified raw materials

Status: 02/06/2024

**EDNA International GmbH**, Gollenhoferstraße 3, 86441 Zusmarshausen, Deutschland,  
Phone: +49 8291 84 390, Email: [export@edna.de](mailto:export@edna.de)